



HIV type 1 vif (IIIB) [FITC] (DAGH008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|----------------------------------------------------------------------------------------------------|
| Antigen Description | E.coli derived vif (HIV type 1) conjugated to FITC. |
| Specificity | This fluoresce conjugated vif binds to murine monoclonal antibodies of defined epitope specificity |
| Conjugate | FITC |
| Applications | Dot Blots; FACS; Immunohistocytology; Direct Immunofluorescence |
| Stability | At least 6 months at -75°C |
| Bio-activity | Not determined |
| Format | Frozen Liquid |
| Size | 50 µg |
| Preservative | None |
| Storage | -75 °C |

BACKGROUND

| | |
|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Keywords | HIV type 1 vif; Viral infectivity factor; HIV-1 vif; HIV1 vif; Human immunodeficiency virus 1; Human immunodeficiency virus type 1 vif; Human immunodeficiency virus 1 Viral infectivity factor |
|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|