



HIV type 1 tat (MN) [FITC] (DAGH012)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Fluorescence conjugated tat HIV type 1 MN (E.coli).
Specificity	FITC-conjugated tat binds to murine monoclonal antibodies of defined epitope specificity and rabbit polyclonal antibodies: dot blot assays detected less than 10ng of tat – long and short UV fluorescent tests.
Conjugate	FITC
Applications	Dot Blot Assays; FITC Cell Staining
Stability	At least 12 months at -75°C
Bio-activity	Not determined.
Size	50 µg
Preservative	None
Storage	-75 °C

BACKGROUND

Keywords	HIV type 1 tat; HIV-1 tat; human immunodeficiency virus-1 tat; human immunodeficiency virus; 30 kDa HIV-1 TAT-interacting protein; HIV-1 TAT-interactive protein 2; Oxidoreductase HTATIP2
-----------------	--