



Aln g 4.0101 (full length) [His] (DAGK-010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Polcalcin Aln g 4 produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 10,185 Dalton. Aln G 4.0101 is expressed with a 6xHis tag at N-terminus and purified by proprietary chromatographic techniques.
Purity	> 80%, based on SDS PAGE
Conjugate	His
Applications	Western blot with human sera.
Molecular Weight	10,185 Dalton
Endotoxin	Tested by LAL (Limulus Amoebocyte Lysate) chromogenic endotoxin assay.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	50 µg, 1 mg
Buffer	Neutral to slightly alkaline pH and 20% glycerol as cryoprotective agent.
Preservative	None
Storage	-70°C or below. Repeated freeze/thaw cycles should be avoided.

BACKGROUND

Introduction	Polcalcin Aln g 4 (Aln G 4.0101) causes an allergic reaction in humans.
Keywords	Polcalcin Aln g 4; Calcium-binding pollen allergen Aln g 4; Aln g 4; Aln g 4.0101