



# Recombinant Human T-Cell Lymphotropic Virus Type II gp21+gp46 chimeric antigen, HRP-Conjugated (DAGRM066T)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	HRP
Applications	Pair for double antigen sandwich assay. For CILa & TRFIA assays
Preservative	None
Storage	Store at +4°C for short term (1-2 weeks) and -20°C for long term. Avoid multiple freeze/thaw cycles.

## BACKGROUND

Introduction	Human T-lymphotropic virus (HTLV) belongs to the family Retroviridae, genus Deltaretrovirus. It is a human, single-stranded RNA retrovirus that causes T-cell leukemia and T-cell lymphoma in adults and may also be involved in certain demyelinating diseases. Adult T-lymphotropic virus (ATLV) is a strain of this disease that effect primarily adults.
Keywords	HTLV-1; gp21; gp46