



Anti-C1QC monoclonal antibody, clone 3R9/2 (CABT-45008MH)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse anti Human C1q antibody, clone 3R9/2 recognizes the human C1q Protein. It is specific for the globular head domain.
Specificity	C1QC
Immunogen	Native C1q
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3R9/2
Conjugate	Unconjugated
Applications	IHC-Fr; ELISA; WB
Format	Liquid
Size	100 µg
Preservative	0.09% Sodium Azide
Storage	-20°C only (ship +4°C)

GENE INFORMATION

Gene Name [C1QC complement component 1, q subcomponent, C chain \[Homo sapiens \(human\) \]](#)

Official Symbol	C1QC
Synonyms	C1QC; complement component 1, q subcomponent, C chain; C1QG; C1Q-C; complement C1q subcomponent subunit C; complement component 1, q subcomponent, gamma polypeptide;
Entrez Gene ID	714
Protein Refseq	NP_001107573
UniProt ID	P02746
Chromosome Location	1p36.11
Pathway	Chagas disease (American trypanosomiasis); Classical antibody-mediated complement activation; Complement Activation, Classical Pathway; Complement and Coagulation Cascades; Complement and coagulation cascades; Complement cascade; Creation of C4 and C2 activators; Immune System;