



Anti-BID monoclonal antibody, clone 25G3 (CABT-B020)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Bid (mouse)
Target	BID
Immunogen	A full length bacterial fusion protein corresponding to Bid
Isotype	IgG2a
Source/Host	Rat
Species Reactivity	Mouse
Clone	25G3
Purification	Protein G affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	1mg/ml
Size	100 µg
Preservative	None
Storage	Store at 4°C
Ship	Cold packs

BACKGROUND

Introduction Bid, a member of BH3 only Bcl-2 family of proteins, transmits the pro-apoptotic signal activated by the death receptor family from the cytoplasmic membrane to the mitochondria. Bid cleaved by caspase-8 to the active form tBid.

Keywords BID;BH3 interacting domain death agonist;FP497;BH3-interacting domain death agonist;p22 BID;BID isoform Si6;BID isoform L(2);BID isoform ES(1b);desmocollin type 4;apoptotic death agonist;Human BID coding sequence

GENE INFORMATION

Entrez Gene ID [637](#)

UniProt ID [A8ASI8](#)