



Mouse anti-Human COL20A1 monoclonal antibody, clone 2G6 (CABT-B10008)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	KIAA1510 (AAH19637, 1 a.a. ~ 156 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2G6
Conjugate	Unconjugated
Applications	sELISA,ELISA
Sequence Similarities	MRGLEGTAGLPGGPRGFQGMAGARGTSGERGPPGTVGPTGLPGPKGERGEKGE PQSLA TLYQLVSQACESAIQTHVSKFDSFHENTRPPMPILEQKLEPGTEPLGSPGTRSKALVPGE WGRGGRHLEGRGEPGAVGQMGSPGQQGASTQGLWE*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	COL20A1 is a collagen gene.Mouse monoclonal antibody raised against a full length
---------------------	---

recombinant KIAA1510.

Keywords	COL20A1; collagen, type XX, alpha 1; bA261N11.4; collagen alpha-1(XX) chain; collagen-like protein;
-----------------	---

GENE INFORMATION

Entrez Gene ID	57642
-----------------------	-----------------------

UniProt ID	Q9P218
-------------------	------------------------

Function	structural molecule activity
-----------------	------------------------------
