



Mouse anti-Human COL24A1 monoclonal antibody, clone 2E9 (CABT-B10010)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	COL24A1 (NP_690850, 1626 a.a. ~ 1715 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2E9
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	PRWTSTQTSGPGLPIGFKGWNGQIFKVNTLLEPKVLSDDCKIQDGSWHKATFLFHTQEPN QLPVEVQKLPPLKTERKYYIDSSVCFL*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords COL24A1; collagen, type XXIV, alpha 1; collagen alpha-1(XXIV) chain;

GENE INFORMATION

Entrez Gene ID	255631
UniProt ID	Q17RW2
Pathway	ATF-2 transcription factor network, organism-specific biosystem
Function	extracellular matrix structural constituent