



Mouse anti-Human COMMD7 monoclonal antibody, clone 2C0 (CABT-B10018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	COMMD7 (AAH47440, 1 a.a. ~ 201 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2C0
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MGRLHCTEDPVPEAVGGDMQQLNQLGAQQFSALTEVLHFHFLTEQKEVERFLAQLSEFATT NQISLGSLRSIVKSLLLVPNGALKKSLTAKQVQADFITLGLSEEKATYFSEKWKQNTPTL ARWAIGQTLMINQLIDMEWKFGVTSGSSELEKVGSI FLQLKLVVKKGNQTKNVYIELTLP QFYSFLHEMERVRTSMECF*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	COMMD7; COMM domain containing 7; C20orf92; dJ1085F17.3; COMM domain-containing protein 7;
-----------------	--

GENE INFORMATION

Entrez Gene ID	149951
-----------------------	------------------------

UniProt ID	Q86VX2
-------------------	------------------------

Function	NF-kappaB binding; protein binding
-----------------	------------------------------------
