



Mouse anti-Human CPXM2 monoclonal antibody, clone 4F5 (CABT-B10030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CPXM2 (NP_937791, 190 a.a. ~ 280 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	4F5
Conjugate	Unconjugated
Applications	ELISA
Sequence Similarities	DLQQWIEVDARRLTRFTGVITQGRNSLWLSDWVTSYKVMVSNDSHTWVTVKNGSGDMIFE GNSEKEIPVLNELPVPVMVARYIRINPQSWF*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CPXM2; carboxypeptidase X (M14 family), member 2; CPX2; UNQ676; inactive carboxypeptidase-like protein X2; 4632435C11Rik; carboxypeptidase Hlo; cytosolic
-----------------	---

GENE INFORMATION

Entrez Gene ID	119587
----------------	------------------------

UniProt ID	Q8N436
------------	------------------------

Function	metallocarboxypeptidase activity; zinc ion binding
----------	--
