



Mouse anti-Human CRELD2 monoclonal antibody, clone 3C7 (CABT-B10031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CRELD2 (NP_077300, 162 a.a. ~ 259 a.a) partial recombinant protein with GST tag.
Isotype	IgM
Source/Host	Mouse
Species Reactivity	Human
Clone	3C7
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	GDGSCRCHMGYQGGLCTDCMDGYFSSLRNETHSICTACDESKTCSGLTNRDCGECEVGW VLDEGACVDVDECAAEPSPCSAAQFCKNANGSYTCEE*
Format	Liquid
Size	200 µl
Buffer	In ascites fluid
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Cysteine-rich with EGF-like domain protein 2 is a protein that in humans is encoded by the CRELD2 gene. Mouse monoclonal antibody raised against a partial recombinant CRELD2.
--------------	--

Keywords CRELD2; cysteine-rich with EGF-like domains 2; cysteine-rich with EGF-like domain protein 2;

GENE INFORMATION

Entrez Gene ID	79174
UniProt ID	Q6UXH1
Function	calcium ion binding
