



Mouse anti-Human CUL9 monoclonal antibody, clone 4G8 (CABT-B10052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PARC (NP_055904, 918 a.a. ~ 1026 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	4G8
Conjugate	Unconjugated
Applications	WB, ELISA
Sequence Similarities	MMLLNRYSEPPGSPERAALETPIIQGDGSPELLIRSLVGGPSAELLLDLERVLCREGSP GGAVRPLLKRLQKETQPFLLLRRLDAPGPNKTLSSVLRVITRLLDF*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	CUL9; cullin 9; PARC; H7AP1; cullin-9; CUL-9; Ubch7-associated protein 1; parkin-like cytoplasmic p53 binding protein; p53-associated parkin-like cytoplasmic protein;
-----------------	--

GENE INFORMATION

Entrez Gene ID	23113
UniProt ID	Q8IWT3
Function	ATP binding; metal ion binding; nucleotide binding; protein binding; ubiquitin protein ligase binding; zinc ion binding