



Mouse anti-Human CYorf15B monoclonal antibody, clone 3G6 (CABT-B10056)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	CYorf15B (NP_115965, 82 a.a. ~ 181 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3G6
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	YKVFQIKLERLEKLYKALQIERNELSEKLGILKGQVSVKVADVDLAVPVTHSCADLDSSN MLNTSSKRAPGVHLEADPKGMNEVKCYSKALSTGSPLGI*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords CYORF15B; chromosome Y open reading frame 15B; lipopolysaccharide-specific response 5-like protein;

GENE INFORMATION

Entrez Gene ID [84663](#)

UniProt ID [Q9BZA4](#)
