



Mouse anti-Human DACT3 monoclonal antibody, clone 3B6 (CABT-B10062)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids 89-202 of human DACT3.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	human
Clone	3B6
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS, pH 7.4 (0.1% sodium azide)
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DACT3; dishevelled-binding antagonist of beta-catenin 3; RRR1; DAPPER3; dapper homolog 3; arginine-rich region 1 protein; antagonist of beta-catenin Dapper homolog 3; dapper, antagonist of beta-catenin, homolog 3;
----------	---

GENE INFORMATION

Entrez Gene ID [147906](#)

UniProt ID [Q96B18](#)
