



Mouse anti-Human DDIT4L monoclonal antibody, clone 3F7 (CABT-B10074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DDIT4L (NP_660287, 94 a.a. ~ 194 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3F7
Conjugate	Unconjugated
Applications	sELISA,ELISA
Sequence Similarities	RLSSTEPCLRGCVMHVNLEIENVCKKLDRIVCDSSVVPTFELTLVFKQENCWSFRDF FFSRGRFSSGFRRTLILSSGFRLVKKKLYSLIGTTVIEGS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DDIT4L; DNA-damage-inducible transcript 4-like; REDD2; Rtp801L; DNA damage-inducible transcript 4-like protein; REDD-2; homolog of mouse SMHS1; HIF-1 responsive protein
-----------------	--

RTP801-like; regulated in development and DNA damage response 2; protein regulated in development and DNA damage response 2;

GENE INFORMATION

Entrez Gene ID [115265](#)

UniProt ID [Q96D03](#)
