



Mouse anti-Human DDX43 monoclonal antibody, clone 4H23 (CABT-B10078)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DDX43 (NP_061135, 2 a.a. ~ 91 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	4H23
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	SHHGGAPKASTWVVASRRSSTVSRAPERRPAEELNRTGPEGYSVGRGGRWRGTSRPPEAV AAGHEELPLCFALKSHFVGAVIGRGSKI*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene is an ATP-dependent RNA helicase in the DEAD-box family and displays tumor-specific expression. [provided by RefSeq, Jul 2008]
---------------------	---

Keywords	DDX43; DEAD (Asp-Glu-Ala-Asp) box polypeptide 43; CT13; HAGE; probable ATP-dependent RNA helicase DDX43; helical antigen; DEAD box protein 43; DEAD-box protein 43; DEAD box protein HAGE; cancer/testis antigen 13;
-----------------	--

GENE INFORMATION

Entrez Gene ID	55510
UniProt ID	Q9NXZ2
Function	ATP binding; ATP-dependent RNA helicase activity; RNA binding; helicase activity; hydrolase activity; nucleotide binding
