



Mouse anti-Human DEPDC1B monoclonal antibody, clone 3I3 (CABT-B10083)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DEPDC1B (NP_060839, 430 a.a. ~ 530 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	3I3
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	SFCRQISPEEFYQRSYGSQEPLAALLEEVITDAKLSNKEKKKKLKQFQKSYPEVYQERF PTPESAALLFPEKPKPKPQLLMWALKKPFQPFQRTSRFRM*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DEPDC1B; DEP domain containing 1B; XTP1; BRCC3; DEP domain-containing protein 1B; breast cancer cell 3; HBxAg transactivated protein 1; HBV XAg-transactivated protein 8; HBV
-----------------	---

GENE INFORMATION

Entrez Gene ID	55789
----------------	-----------------------

UniProt ID	Q8WUY9
------------	------------------------

Pathway	Rho GTPase cycle, organism-specific biosystem; Signaling by Rho GTPases, organism-specific biosystem
---------	--

Function	GTPase activator activity
----------	---------------------------
