



Mouse anti-Human DLX6 monoclonal antibody, clone 3E8 (CABT-B10107)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DLX6 (XP_376652, 71 a.a. ~ 161 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3E8
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	PYASGGGNSYNHRSLAAYPYMSHSQHSPLYLSYHNSSAAAQTRGDDTDQKKTTVIENGEI RFNGKGKKIRKPRTIYSSLQLQALNHRFQQ*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene encodes a member of a homeobox transcription factor gene family similar to the Drosophila distal-less gene. This family is comprised of at least 6 different members that
---------------------	---

encode proteins with roles in forebrain and craniofacial development. This gene is in a tail-to-tail configuration with another member of the family on the long arm of chromosome 7.
[provided by RefSeq, Jul 2008]

Keywords	DLX6; distal-less homeobox 6; homeobox protein DLX-6; distal-less homeo box 6;
-----------------	--

GENE INFORMATION

Entrez Gene ID	1750
-----------------------	----------------------

UniProt ID	P56179
-------------------	------------------------

Function	sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity
-----------------	--
