



Mouse anti-Human DMWD monoclonal antibody, clone 4G6 (CABT-B10113)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DMWD (AAH19266, 245 a.a. ~ 334 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	4G6
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	GGEKPSGPVPRSRLDPAKVLGTALCPRIHEVPLLEPLVCKKIAQERLTVLLFLEDCIITA CQEGLECTWARPGKAGISSQPGNSPSGTVV
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Dystrophia myotonica WD repeat-containing protein is a protein that in humans is encoded by the DMWD gene. Mouse monoclonal antibody raised against a partial recombinant DMWD.
---------------------	---

Keywords

DMWD; dystrophia myotonica, WD repeat containing; DMRN9; DMR-N9; gene59; D19S593E; dystrophia myotonica WD repeat-containing protein; protein 59; dystrophia myotonica-containing WD repeat motif protein;

GENE INFORMATION

Entrez Gene ID	1762
UniProt ID	Q8WUW6
Function	molecular_function