



# Mouse anti-Human DPY30 monoclonal antibody (CABT-B10133)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to full length human DPY30.
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB,ELISA
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS, pH 7.4 (0.02% sodium azide, 1% BSA, 50% glycerol)
<b>Storage</b>	At -20°C for one year.

## BACKGROUND

<b>Introduction</b>	DPY30 is a component of a Set1-like multiprotein histone methyltransferase complex (Cho et al., 2007 [PubMed 17500065]).[supplied by OMIM, May 2008]
<b>Keywords</b>	DPY30; dpy-30 homolog (C. elegans); Cps25; Saf19; HDPY-30; protein dpy-30 homolog; dpy-30L; dpy-30-like protein;

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">84661</a>
<b>UniProt ID</b>	<a href="#">Q9C005</a>
<b>Function</b>	identical protein binding; protein binding; protein homodimerization activity