



Mouse anti-Human DZIP3 monoclonal antibody, clone 4B2 (CABT-B10156)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	DZIP3 (NP_055463, 1021 a.a. ~ 1121 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	4B2
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	PSAGLRSDPSIMNWERITDRLKTAFPQQTRKELTDFLRKLKDAYGKSLSELTfDEIVCKI SQFIDPKKSQSQGKSVSNVNCVSPSHSPSQPDAAQPPKPA*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	DZIP3; DAZ interacting zinc finger protein 3; UURF2; PPP1R66; hRUL138; E3 ubiquitin-protein ligase DZIP3; UURF2 ubiquitin ligase; DAZ-interacting protein 3; zinc finger DAZ interacting
-----------------	--

protein 3; DAZ interacting protein 3, zinc finger; RNA-binding ubiquitin ligase of 138 kDa;
protein phosphatase 1, regulatory subunit 66; human RNA-binding ubiquitin ligase of 138 kDa;

GENE INFORMATION

Entrez Gene ID	9666
UniProt ID	Q86Y13
Pathway	Adaptive Immunity Signaling, organism-specific biosystem; Antigen processing: Ubiquitination & Proteasome degradation, organism-specific biosystem; Class I MHC mediated antigen processing & presentation, organism-specific biosystem; Immune System, organism-specific biosystem
Function	RNA binding; ligase activity; metal ion binding; polyubiquitin binding; ubiquitin-protein ligase activity; zinc ion binding