



Mouse anti-Human EIF2C4 monoclonal antibody, clone 5H4 (CABT-B10175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	EIF2C4 (NP_060099, 1 a.a. ~ 101 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5H4
Conjugate	Unconjugated
Applications	sELISA, ELISA
Sequence Similarities	MEALGPGPPASLFQPPRRPGLGTVGKPIRLLANHFQVQIPKIDVYHYDVDIKPEKRPRRV NREVVDTMVRHFKMQIFGDRQPGYDGKRNMYYAHPLPIGR*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	eukaryotic translation initiation factor 2C, 4; FLJ20033; AGO4; argonaute 4; KIAA1567; argonaute4; protein argonaute-4; hAgo4; Eukaryotic translation initiation factor 2C 4; eIF2C 4;
-----------------	--

GENE INFORMATION

Entrez Gene ID	192670
----------------	------------------------

UniProt ID	Q9HCK5
------------	------------------------

Pathway	MicroRNA (miRNA) Biogenesis, organism-specific biosystem; Regulatory RNA pathways, organism-specific biosystem
---------	--

Function	RNA binding; protein binding
----------	------------------------------
