



Mouse anti-Mouse env monoclonal antibody, clone 25-5-5T (CABT-B10194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	An epitope on mouse I-Ek MHC class II alloantigen.
Immunogen	Native from C3H mouse skin graft and splenocytes.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Mouse
Clone	25-5-5T
Conjugate	Unconjugated
Applications	IHC, IP, FC
Format	Liquid
Buffer	In 100 mM BBS, pH 8.2
Storage	Store at 4°C.

BACKGROUND

Keywords ENV; gp160; envelope glycoprotein; Envelope surface glycoprotein gp160, precursor; hypothetical protein; Envelope surface glycoprotein gp120; Envelope transmembrane domain;

GENE INFORMATION

