



Mouse anti-Human ETV5 monoclonal antibody, clone 4C21 (CABT-B10210)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ETV5 (NP_004445, 181 a.a. ~ 291 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	4C21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	LPEPGPQQQTFAVPRPPHQPLQMPKMMMPENQYPSEQRFRQLSEPCHPFPQPGVPGDNR PSYHRQMSEPIVPAAPPPPQGFKQEYHDPLEYHGVPGMPGPPAHGFQSPM*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	ETV5; ets variant 5; ERM; ETS translocation variant 5; ets-related molecule; ets-related protein ERM;
-----------------	---

GENE INFORMATION

Entrez Gene ID	2119
UniProt ID	P41161
Pathway	Androgen Receptor Signaling Pathway, organism-specific biosystem; FOXM1 transcription factor network, organism-specific biosystem; IL12 signaling mediated by STAT4, organism-specific biosystem
Function	DNA binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; transcription regulatory region DNA binding