



Mouse anti-Human FAIM2 monoclonal antibody, clone 3D3 (CABT-B10223)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FAIM2 (AAH00051, 1 a.a. ~ 317 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3D3
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MTQGKLSVANKAPGTEGQQQVHGEKKEAPAVPSAPPSYEEATSGEGMKAGAFPPAPTAVP LHPSWAYVDPSSSSSYDNGFPTGDHELFTTFSWDDQKVRRVFVRKVYTILLIQLLVTLAV VALFTFCDPVKDYVQANPGWYWASYAVFFATYLTACCSPRRHFPWNLILLTVFTLSMA YLTGMLSSYYNTTSVLLCLGITALVCLSVTVFSFQTKFDFTSCQGVLFVLLMTLFFSGLI LAILLPFQYVPWLHA
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	FAIM2; Fas apoptotic inhibitory molecule 2; LFG; LFG2; NGP35; NMP35; TMBIM2; protein lifeguard 2; neural membrane protein 35; transmembrane BAX inhibitor motif-containing protein 2;
-----------------	---

GENE INFORMATION

Entrez Gene ID	23017
-----------------------	-----------------------

UniProt ID	Q9BWQ8
-------------------	------------------------

Pathway	FAS (CD95) signaling pathway, organism-specific biosystem
----------------	---
