



Mouse anti-Human FAIM3 monoclonal antibody, clone 2F5 (CABT-B10224)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FAIM3 (AAH06401, 124 a.a. ~ 223 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2F5
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	SEYEPSWEEQMPETPKWFHLPYLFQMPAYASSSKFVTRVTTPAQRGKVPPVHHSSPTTQ ITHRPRVSRASSVAGDKPRTFLPSTTASKISALEGLLKPQ
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Fc receptors specifically bind to the Fc region of immunoglobulins (Igs) to mediate the unique functions of each Ig class. FAIM3 encodes an Fc receptor for IgM (see MIM 147020)
---------------------	--

(Kubagawa et al., 2009 [PubMed 19858324]; Shima et al., 2010 [PubMed 20042454]).[supplied by OMIM, Jul 2010]

Keywords

FAIM3; Fas apoptotic inhibitory molecule 3; FCMR; TOSO; fas apoptotic inhibitory molecule 3; Fc mu receptor; IgM Fc receptor; immunoglobulin mu Fc receptor; regulator of Fas-induced apoptosis Toso;

GENE INFORMATION

Entrez Gene ID

[9214](#)

UniProt ID

[O60667](#)
