



Mouse anti-Human FAM65B monoclonal antibody, clone 3F72C21 (CABT-B10226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	C6orf32 (AAH01232, 1 a.a. ~ 424 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	3F72C21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MLVGSQSFSPPGGPNGIIRSQSFAGFSGLQERRSRCNSFIENSSALKKPQAKLKKMHNLGH KNNNPPKEPQPKRVEEVYRALKNGLDEYLEVHQTELDKLTATQLKDMKRNSRLGVLYDLDK QIKTIERYMRRLEFHISKVDELYEAYCIQRRQLQDGASKMKQAFATSPASKAARESLTEIN RSFKEYTENMCTIEVELENLLGEFSIKMKGLAGFARLCPGDQYEIFMKYGRQRWKLKGKI EVNGKQSWDGEETVF
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene stimulates the formation of a non-mitotic multinucleate syncytium from proliferative cytotrophoblasts during trophoblast differentiation. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Nov 2013]
Keywords	FAM65B; family with sequence similarity 65, member B; PL48; DIFF40; DIFF48; MYONAP; C6orf32; protein FAM65B; myogenesis-related and NCAM-associated protein homolog;

GENE INFORMATION

Entrez Gene ID	9750
UniProt ID	Q9Y4F9
Function	binding