



Mouse anti-Human FAM84A monoclonal antibody, clone 2D3 (CABT-B10227)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NSE1 (AAH26346, 1 a.a. ~ 293 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2D3
Conjugate	Unconjugated
Applications	WB,IHC,sELISA,ELISA
Sequence Similarities	MGNQLDRITHLNYSELPTGDPGIEKDELRVGVAYFFSDDEEDLDERGQPDKFGVKAPPG CTPCPESPSRHHHLLHQLVLNETQFSAFRGQECIFSKVSGGPQGADLSVYAVTALPALC EPGDLELLWLQHAPEPPAPAPHWAVYVGGGQIIHLHQGEIRQDSL YEAGAANVGRVVNS WYRYRPLVAELVVQNACGHLGLKSEEICWTNSESFAAWCRFGKREFKAGGEVPAGTQPPQ QQYYLKVHLGENKVH
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords FAM84A; family with sequence similarity 84, member A; NSE1; PP11517; protein FAM84A; neurological/sensory 1; neurologic sensory protein 1;

GENE INFORMATION

Entrez Gene ID [151354](#)

UniProt ID [Q96KN4](#)
