



Mouse anti-Human FCHSD1 monoclonal antibody, clone 5H5 (CABT-B10254)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FCHSD1 (AAH47016, 1 a.a. ~ 410 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5H5
Conjugate	Unconjugated
Applications	ELISA
Sequence Similarities	MDSRTVFGAWRCLLDATVAGGQTRLQASDRYRDLAGGTGRSAKEQVLRKGTENLQRAQAE VLQSVRELSRSRKLYGQRERWALAEKAADVQARLNRSDHGIFHSRTSLQKLSTKLSAQ SAQYSQQLQAARNEYLLNLVATNAHLHDHYYQEELPALLKALVSELSEHLRDPLTSLSHT LEAAEVILEHAHRGEQTTSQVSWEQDLKLFLQEPGVFSPTPPQQFQFQAGTDQVCVLEWGA EGVAGKSGLEKEVQR
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	FCH and double SH3 domains protein 1 is a protein that in humans is encoded by the FCHSD1 gene. Mouse monoclonal antibody raised against a full length recombinant FCHSD1.
Keywords	FCHSD1; FCH and double SH3 domains 1; NWK2; FCH and double SH3 domains protein 1; nervous wreck homolog 2;

GENE INFORMATION

Entrez Gene ID	89848
UniProt ID	Q86WN1