



Mouse anti-Human FGFR1OP2 monoclonal antibody, clone 3H5 (CABT-B10264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FGFR1OP2 (NP_056448, 62 a.a. ~ 170 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3H5
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	RSTLVMGIQQENRQIRELQQENKELRTSLEEHSALIMSKYREQMFRLLMASKKDDPG IIMKLKEQHISKIDMVHRNKSEGFFLDASRHILEAPQHGLERRHLEANQ*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	FGFR1OP2; FGFR1 oncogene partner 2; WIT3.0; HSPC123-like; wound inducible transcript 3.0; fibroblast growth factor receptor 1 oncogene partner 2;
-----------------	---

GENE INFORMATION

Entrez Gene ID [26127](#)

UniProt ID [Q9NVK5](#)
