



Mouse anti-Human FMN2 monoclonal antibody, clone 5C9 (CABT-B10276)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FMN2 (AAH14364, 144 a.a. ~ 243 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	5C9
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	*PRGRGKFTDRDS*MLFGDHGIFLHETKTWREGGVPKCFLQYLA*IQL*L*RLLEEREQT SSTRESKRSRRGV*TEERKITL*NKTKT*LWN*SKDKHEN
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Formin homology (FH) domain proteins (see FMN1; MIM 136535) play a role in cytoskeletal organization and/or establishment of cell polarity.[supplied by OMIM, Apr 2010]
---------------------	---

Keywords

FMN2; formin 2; formin-2;

GENE INFORMATION

Entrez Gene ID[56776](#)

UniProt ID[Q9NZ56](#)

Pathway

Dorso-ventral axis formation, organism-specific biosystem; Dorso-ventral axis formation, conserved biosystem

Function

actin binding; molecular_function