



Mouse anti-Human FUT10 monoclonal antibody, clone 5I4 (CABT-B10287)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FUT10 (AAH04884, 1 a.a. ~ 93 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	5I4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MEEAPTHLNSFLKKEGLTFNRKRKWELDSYPIMLWWSPLTGETGRLGQCGADACFFTINR TYLHHHMTKAFLFYGKQDFRLSPLFAVFLQS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	FUT10; fucosyltransferase 10 (alpha (1,3) fucosyltransferase); alpha-(1,3)-fucosyltransferase 10; fuc-TX; fucT-X; fucosyltransferase X; alpha 1,3-fucosyl transferase; galactoside 3-L-
-----------------	---

GENE INFORMATION

Entrez Gene ID [84750](#)

UniProt ID [Q6P4F1](#)

Function alpha(1,3)-fucosyltransferase activity; transferase activity, transferring glycosyl groups
