



Mouse anti-Human FXC1 monoclonal antibody, clone 2B22 (CABT-B10290)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FXC1 (NP_036324, 11 a.a. ~ 104 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2B22
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	LRNLRDFLLVYNRMTELCFQRCVPSLHHRALDAEEEEACLHSCAGKLIHSNHRLMAAYVQL MPALVQRRRIADYEAASAVPGVAAEQPGVSPSGS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	fracture callus 1 homolog (rat); OTTHUMP00000213157; TIMM10B; mitochondrial import inner membrane translocase subunit Tim9 B; TIM10B; Tim9b; Fracture callus protein 1; TIMM9B;
-----------------	---

GENE INFORMATION

Entrez Gene ID	26515
----------------	-----------------------

UniProt ID	Q9Y5J6
------------	------------------------

Function	metal ion binding
----------	-------------------
