



# Rat anti-Mouse Gen1 monoclonal antibody, clone USB09 (CABT-B10318)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	The hybridoma producing USB09 was obtained by immunization of a rat with the cell lysate of adult mouse testis.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rat
<b>Species Reactivity</b>	Mouse
<b>Clone</b>	USB09
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In PBS (50% glycerol)
<b>Storage</b>	Store at -20°C. For long term storage store at -80°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Keywords</b>	GEN1; Gen homolog 1, endonuclease (Drosophila); 5830483C08Rik; flap endonuclease GEN homolog 1;
-----------------	---

## GENE INFORMATION

**Entrez Gene ID** [209334](#)

---

**Function** DNA binding; catalytic activity; endonuclease activity; hydrolase activity; metal ion binding; nuclease activity;

---