



Mouse anti-Human GK5 monoclonal antibody, clone 3D22 (CABT-B10330)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MGC40579 (NP_689989, 1 a.a. ~ 90 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3D22
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	MSGLLTDPEQRAQEPRYPGFVLGLDVGSSVIRCHVYDRAARVCGSSVQKVENLYPQIGWV EIDPDVLWIQFVAVIKEAVKAAGIQMNQI*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords GK5; glycerol kinase 5 (putative); putative glycerol kinase 5; GK 5; glycerokinase 5; ATP:glycerol 3-phosphotransferase 5;

GENE INFORMATION

Entrez Gene ID [256356](#)

UniProt ID [Q6ZS86](#)

Function ATP binding; glycerol kinase activity; nucleotide binding; transferase activity
