



Mouse anti-Human GLTSCR2 monoclonal antibody, clone 6B9 (CABT-B10337)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GLTSCR2 (NP_056525, 402 a.a. ~ 479 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	6B9
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	KPRRLGRLKYQAPDIDVQLSSELTDSLRTLKPEGNILRDRFKSFQRRNMIEPRERAKFKR KYKVKLVEKRAFREIQL*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Glioma tumor suppressor candidate region gene 2 protein is a protein that in humans is encoded by the GLTSCR2 gene. Mouse monoclonal antibody raised against a partial
---------------------	--

recombinant GLTSCR2.

Keywords

GLTSCR2; glioma tumor suppressor candidate region gene 2; PICT1; PICT-1; glioma tumor suppressor candidate region gene 2 protein; p60; protein interacting with carboxyl terminus 1;

GENE INFORMATION

Entrez Gene ID

[29997](#)

UniProt ID

[Q9NZM5](#)

Function

molecular_function
