



Mouse anti-Human GMFB monoclonal antibody, clone 3H233B3 (CABT-B10338)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GMFB (AAH05359, 1 a.a. ~ 155 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3H233B3
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MSESLVVCDAEDLVEKLRKFRFRKETNNAIIMKIDKDKRLVVLDEELEGISPDELKDE LPERQPRFIVYSYKYQHDDGRVSYPLCFIFSSPVGCKPEQQMMYAGSKNKLVTAEATKV FEIRNTEDLTEEWLREKLGFFTNVNFVSKVFMV*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	GMFB; glia maturation factor, beta; GMF; glia maturation factor beta; GMF-beta;
-----------------	---

GENE INFORMATION

Entrez Gene ID	2764
UniProt ID	P60983
Function	actin binding; enzyme activator activity; growth factor activity; protein kinase inhibitor activity; signal transducer activity