



Mouse anti-Human GOLGA7 monoclonal antibody, clone 3I9 (CABT-B10351)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GOLGA7 (AAH12032, 1 a.a. ~ 138 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	3I9
Conjugate	Unconjugated
Applications	WB,IP,sELISA,ELISA
Sequence Similarities	MRPQQAPVSGKVFIQRDYSSGTRCQFQTKFPAELENRIDRQQFEETVRTLNNLYAEAEKL GGQSYLEGCLACLTAYTIFLCMETHYEKVLKKVSKYIQEQNEKIYAPQGLLLTDPPIERGL RVIEITIYEDRGMSSGR*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords GOLGA7; golgin A7; GCP16; GOLGA7A; HSPC041; GOLGA3AP1; Golgin subfamily A

member 7; golgi autoantigen, golgin subfamily a, 7; Golgi complex-associated protein of 16kDa;
golgi complex-associated protein of 16 kDa;

GENE INFORMATION

Entrez Gene ID [51125](#)

UniProt ID [Q7Z5G4](#)
