



Mouse anti-Human GOLGB1 monoclonal antibody, clone 7E5 (CABT-B10352)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GOLGB1 (NP_004478, 3 a.a. ~ 92 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	7E5
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	SRLSGLANVVLHELSGDDDDTDQNMRAPLDPELHQESDMEFNNTTQEDVQERLAYAEQLVV ELKDIIRQKDVQLQKDEALQEERKAADN*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Golgin subfamily B member 1 is a protein that in humans is encoded by the GOLGB1 gene. Mouse monoclonal antibody raised against a partial recombinant GOLGB1.
---------------------	---

Keywords	GOLGB1; golgin B1; GCP; GCP372; GOLIM1; Golgin subfamily B member 1; giantin; macrogolgin; golgi integral membrane protein 1; 372 kDa Golgi complex-associated protein; golgin B1, golgi integral membrane protein; golgi autoantigen, golgin subfamily b, macrogolgin (with transmembrane signal), 1;
-----------------	--

GENE INFORMATION

Entrez Gene ID	2804
UniProt ID	Q14789
Function	protein binding
