



Mouse anti-Human GPR84 monoclonal antibody, clone 6D4 (CABT-B10356)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	GPR84 (NP_065103, 208 a.a. ~ 317 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	6D4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	AAQALDQYKLRQASIHSNHVARTDEAMPGRFQELDSRLASGGPSEGISEPVSAATTQTL EGDSSEVGDQINSKRAKQMAEKSPPEASAKAQPIKGARRAPDSSEFGK*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	GPR84; G protein-coupled receptor 84; EX33; GPCR4; G-protein coupled receptor 84; inflammation-related G protein-coupled receptor EX33; inflammation-related G-protein coupled
-----------------	--

receptor EX33;

GENE INFORMATION

Entrez Gene ID	53831
----------------	-----------------------

UniProt ID	Q9NQS5
------------	------------------------

Pathway	GPCRs, Other, organism-specific biosystem
---------	---

Function	G-protein coupled receptor activity; receptor activity
----------	--
