



Mouse anti-Human HBXIP monoclonal antibody, clone 5H2 (CABT-B10387)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HBXIP (AAH62619, 83 a.a. ~ 174 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5H2
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MEATLEQHLEDTMKNPSIVGVLC TDSQGLNLGCRGTL SDEHAGVISVLAQQA AKLTSDPT DIPVVCLESDNGNIMI QKHDGITVAAHKMAS*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	hepatitis B virus x interacting protein; HBV X-interacting protein; XIP; hepatitis B virus X-interacting protein; MGC71071; hepatitis B virus x-interacting protein (9.6kD); HBX-interacting
-----------------	--

GENE INFORMATION

Entrez Gene ID	10542
----------------	-----------------------

UniProt ID	O43504
------------	------------------------

Function	protein binding
----------	-----------------
