



Mouse anti-Human HECTD1 monoclonal antibody, clone 2F21 (CABT-B10392)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HECTD1 (NP_056197, 3 a.a. ~ 111 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	2F21
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	DVDPDTLLEWLQMGQGDERDMQLIALEQLCMLLLMSDNVDRCFETCPPRTFLPALCKIFL DESAPDNVLEVTARAITYYLDVSAECTRRIVGVGDGAIKALCNRLVVVE*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	HECTD1; HECT domain containing E3 ubiquitin protein ligase 1; EULIR; E3 ubiquitin-protein ligase HECTD1; E3 ligase for inhibin receptor; HECT domain-containing protein 1;
-----------------	--

GENE INFORMATION

Entrez Gene ID	25831
UniProt ID	Q9ULT8
Function	acid-amino acid ligase activity; ligase activity; metal ion binding; protein binding; ubiquitin-protein ligase activity