



Mouse anti-Human HES2 monoclonal antibody, clone 5I7 (CABT-B10396)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HES2 (NP_061962, 41 a.a. ~ 114 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	5I7
Conjugate	Unconjugated
Applications	WB,IF,ELISA
Sequence Similarities	LPLLGRENSNCSKLEKADVLEMTVRFLQELPASSWPTAAPLPDCSYREGYSACVARLARV LPACRVLEPAVSAR
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	HES2; hes family bHLH transcription factor 2; bHLHb40; transcription factor HES-2; hairy and enhancer of split 2; class B basic helix-loop-helix protein 40;
-----------------	--

GENE INFORMATION

Entrez Gene ID [54626](#)

UniProt ID [Q9Y543](#)

Function DNA binding; double-stranded DNA binding
