



Mouse anti-Human HHLA3 monoclonal antibody, clone 2G7 (CABT-B10399)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HHLA3 (NP_009002, 1 a.a. ~ 54 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	2G7
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MFGACYKQPLKPSGSEPPAEECRMTPRHAGCDVTEMQRILSQPTFTEHLLRAV*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	HHLA3; HERV-H LTR-associating 3; HERV-H LTR-associating protein 3; human endogenous retrovirus-H long terminal repeat-associating 3; human endogenous retrovirus-H long terminal repeat-associating protein 3;
-----------------	--

GENE INFORMATION

Entrez Gene ID	11147
UniProt ID	Q9XRX5
Function	protein binding