



Mouse anti-Human HS1BP3 monoclonal antibody, clone 2F3 (CABT-B10428)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HS1BP3 (AAH27947, 1 a.a. ~ 214 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2F3
Conjugate	Unconjugated
Applications	WB,ELISA
Sequence Similarities	MQSPAVLVTSRRLQNAHTGLDLTPQHQEVRGKMMSGHVEYQILVVTRLAAFKSAKHRPE DVVQFLVSKKYSEIEEFYQKLSSRYAAASLPPLPRKVLVFGESDIRERRAVFNEILRCVS KDAELAGSPELLEFLGTRSPGAAGLTSRDSSVLDGTDSTGNDEEAFDFFEEQDQVAEEG PPVQSLKGEDAEESELEEEALDPLGIMRLVSCC*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation. [provided by RefSeq, Jul 2008]
Keywords	HS1BP3; HCLS1 binding protein 3; ETM2; HS1-BP3; HCLS1-binding protein 3; HSP1BP-3; HS1-binding protein 3;

GENE INFORMATION

Entrez Gene ID	64342
UniProt ID	Q9H7U7
Function	phosphatidylinositol binding