



Mouse anti-Human HUNK monoclonal antibody, clone 3F4 (CABT-B10438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	HUNK (NP_055401, 615 a.a. ~ 714 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b
Source/Host	Mouse
Species Reactivity	Human
Clone	3F4
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	LHPTLVSF AHEDKNSPPKEEGLCCPPPVSNGPMQPLGSPNCVKSRGRFPMMGIGQMLRK RHQSLQPSADRPLEASLPPLQPLAPVNLAFDMADGVKTQC
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	HUNK; hormonally upregulated Neu-associated kinase; Bstk1; Mak-v; hormonally up-regulated neu tumor-associated kinase; serine/threonine-protein kinase MAK-V;
-----------------	---

GENE INFORMATION

Entrez Gene ID [30811](#)

UniProt ID [P57058](#)

Function ATP binding; nucleotide binding; protein serine/threonine kinase activity
