



Mouse anti-Human IFRD2 monoclonal antibody, clone 2B52H2 (CABT-B10448)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	IFRD2 (AAH01327, 1 a.a. ~ 443 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	2B52H2
Conjugate	Unconjugated
Applications	WB,IF,sELISA,ELISA
Sequence Similarities	MPRARKGNTLRKGGQRRGGGARSSAQADSGSSDDEAASEARSTASECPSLLSTTAEDSLG GDVVDEQGQQEDLEEKLKEYVDCLTDKSAKTRQGALESLRLALASRLLPDFLLERRRLTLA DALEKCLKKGKGEEQALAAAVLGLLCVQLGPGPKGEELFHSQPLLVSFLSDSTASPAAR LHCASALGLGCYVAAADIQDLVSCLACLESVFSRFYGLGGSSTSPVVPASLHGLLSAALQ AWALLLTICPSTQIS
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	IFRD2; interferon-related developmental regulator 2; SM15; IFNRP; SKMc15; interferon-related protein;
-----------------	---

GENE INFORMATION

Entrez Gene ID	7866
-----------------------	----------------------

UniProt ID	Q12894
-------------------	------------------------

Function	binding; molecular_function
-----------------	-----------------------------
