



Mouse anti-Human IGDCC4 monoclonal antibody, clone MAAG1058 (CABT-B10450)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NOPE (NP_066013, 1152 a.a. ~ 1251 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a
Source/Host	Mouse
Species Reactivity	Human
Clone	MAAG1058
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Clear, colorless solution in phosphate buffered saline, pH 7.2
Size	100 µg
Buffer	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

BACKGROUND

Keywords	IGDCC4; immunoglobulin superfamily, DCC subclass, member 4; Nope; DDM36; immunoglobulin superfamily DCC subclass member 4; amplimer WI-16786, mouse homolog; amplimer WI-18508, mouse homolog; neighbor of Punc e11 protein;
-----------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

GENE INFORMATION

Entrez Gene ID [57722](#)

UniProt ID [Q8TDY8](#)
