



Mouse anti-Human IGSF6 monoclonal antibody, clone 3B23 (CABT-B10452)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	IGSF6 (AAH17844, 29 a.a. ~ 242 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3B23
Conjugate	Unconjugated
Applications	WB,sELISA,ELISA
Sequence Similarities	TLSVTQPWYLEVDYTHEAVTIKCTFSATGCPSEQPTCLWFRYGAHRPENLCLDGCKSEAD KFTVREALKENQVSLTVNRVTSNDSAIYICGIAFPSVPEARAKQTGGGTTLVVREIKLLS KELRNFLTALVSLLSVYVTGVCVAFILLSKSKSNPLRNKKIKEDSQKKKSARRIFQEIAQ ELYHKRHVETNQQSEKDNNTYENRRVLSNYERP*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	IGSF6; immunoglobulin superfamily, member 6; DORA; immunoglobulin superfamily member 6; down-regulated by activation (immunoglobulin superfamily);
-----------------	--

GENE INFORMATION

Entrez Gene ID	10261
-----------------------	-----------------------

UniProt ID	O95976
-------------------	------------------------

Function	transmembrane receptor activity
-----------------	---------------------------------
