



Mouse anti-Human JMY monoclonal antibody, clone 3E8 (CABT-B10483)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	JMY (NP_689618, 535 a.a. ~ 635 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	3E8
Conjugate	Unconjugated
Applications	WB,IF,ELISA
Sequence Similarities	SPMDEVLASLKRGSFHLKKVEQRTLPPFPDEDDSNILAQIRKGVKLKKVQKDVLRRESFT LLPDTDPLTRSIHEALRRRIKEASPESEDEEEEALPCTDWEN*
Format	Liquid
Size	100 µg
Buffer	In 1x PBS, pH 7.2
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	JMY; junction-mediating and regulatory protein; junction-mediating and -regulatory protein; p300 transcriptional cofactor JMY;
-----------------	--

GENE INFORMATION

Entrez Gene ID	133746
UniProt ID	Q8N9B5
Pathway	Direct p53 effectors, organism-specific biosystem
Function	actin binding; protein binding; transcription coactivator activity